



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Cold box,BK

**Product Code :**  
LBNY-0001-2100020



**Description :**

Cold box, Small, Long range  
Vaccine storage capacity: 12 litres  
Cold life at 43°C without openings: 156 hours  
With hinged lid and carrying handles.  
Supplied with 42 x 0.4L icepacks as a standard.  
Weight fully loaded: 45 Kg

**Technical Specification :**

Dimensions:  
Internal dimensions (LxWxH cm): 34 x 33 x 34.8  
External dimensions (LxWxH cm): 61 x 60 x 56



**Laboratory Instrument India**