



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cold box,BK

Product Code :
LBNY-0001-2100020



Description :

Cold box, Small, Long range
Vaccine storage capacity: 12 litres
Cold life at 43°C without openings: 156 hours
With hinged lid and carrying handles.
Supplied with 42 x 0.4L icepacks as a standard.
Weight fully loaded: 45 Kg

Technical Specification :

Dimensions:
Internal dimensions (LxWxH cm): 34 x 33 x 34.8
External dimensions (LxWxH cm): 61 x 60 x 56



Laboratory Instrument India