



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Vacc. carrier,BK

**Product Code :**  
LBNY-0001-2100021



**Description :**

Vaccine carrier, Extra-large, Long range  
Vaccine storage capacity: 3.4 litres  
With lid and shoulder strap.  
Cold life at 43°C without openings: 43.43 hours  
Supplied with 4 x of 0.6L icepacks as a standard.  
Weight fully loaded: 4.9 Kg.

**Technical Specification :**

Dimensions:  
Internal dimensions (LxWxH cm): 13.3 x 1 3.3 x 19  
External dimensions (LxWxH cm): 28.8 x 28.9 x 33.7



**Laboratory Instrument India**