



Email : sales@laboratoryinstrumentindia.com

Product Name :
Vacc. carrier,BK

Product Code :
LBNY-0001-2100021



Description :

Vaccine carrier, Extra-large, Long range
Vaccine storage capacity: 3.4 litres
With lid and shoulder strap.
Cold life at 43°C without openings: 43.43 hours
Supplied with 4 x of 0.6L icepacks as a standard.
Weight fully loaded: 4.9 Kg.

Technical Specification :

Dimensions:
Internal dimensions (LxWxH cm): 13.3 x 1 3.3 x 19
External dimensions (LxWxH cm): 28.8 x 28.9 x 33.7



Laboratory Instrument India