



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Vacc. carrier 0.9 litres

**Product Code :**  
LBNY-0001-210022



**Description :**

Vaccine carrier, Small, Short range  
Vaccine storage capacity: 0.9 litres  
With lid and shoulder strap.  
Cold life at 43°C without openings: 20 hours  
Supplied with 2 x 0.4L icepacks as a standard.  
Weight fully loaded: 2.2 Kg

**Technical Specification :**

Dimensions:  
Internal dimensions (LxWxH cm): 13 x 10 x 16.5  
External dimensions (LxWxH cm): 22.5 x 16 x 25



**Laboratory Instrument India**