



Email : sales@laboratoryinstrumentindia.com

Product Name :
Vacc. carr,Nilkamal

Product Code :
LBNY-0001-2100023



Description :

Vaccine carrier, Extra Large, Long range
Vaccine storage capacity: 2.7 litres
With lid and shoulder strap.
Cold life at 43°C without openings: 46 hours 23 minutes
Supplied with 4 x 0.6L icepacks as a standard.
Weight fully loaded: 6.4 Kg

Technical Specification :

Dimensions:
Internal dimensions (LxWxH cm): 19.0 x 19.0 x 19.0
External dimensions (LxWxH cm): 27.0 x 27.0 x 32.0



Laboratory Instrument India