



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cold box, Nilkamal 6 litres

Product Code :
LBNY-0001-2100024



Description :

Cold box, Small, Long range
Vaccine storage capacity: 6 litres
With hinged lid and carrying handles.
Cold life at 43°C without openings: 106 hours
Supplied with 26 x 0.4 litre icepacks as a standard.
Weight fully loaded: 27.54 Kg.

Technical Specification :

Dimensions:
Internal dimensions (LxWxH cm): 40 x 27.8 x 20
External dimensions (LxWxH cm): 65 x 53 x 46



Laboratory Instrument India