



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cold box, Nilkamal 20.25 litres

Product Code :
LBNY-0001-2100025



Description :

Cold box, Large, Long range
Vaccine storage capacity: 20.25 litres
With hinged lid and carrying handles.
Cold life at 43°C without openings: 152 hours 28 minutes
Supplied with 44 x 0.4L icepacks as a standard.
Weight fully loaded: 49.70 Kg

Technical Specification :

Dimensions:
Internal dimensions (LxWxH cm): 51.00 x 36.00 x 27.50
External dimensions (LxWxH cm): 77.00 x 62.00 x 53.50



Laboratory Instrument India