



Email : sales@laboratoryinstrumentindia.com

Product Name :
Vacc. carr, 1.67 litres

Product Code :
LBNY-0001-2100026



Description :

Vaccine carrier, Large, Long range
Vaccine storage capacity: 1.67 litres
With lid and shoulder strap.
Cold life at 43°C without openings: 35.7 hours
Supplied with 4 x 0.4L icepacks as a standard.
Weight fully loaded: 4.34 Kg.

Technical Specification :

Dimensions:
Internal dimensions (LxWxH cm): 10.45 x 10.45 x 1.67
External dimensions (LxWxH cm): 25 x 25 x 30



Laboratory Instrument India