



Email : sales@laboratoryinstrumentindia.com

Product Name :
Vaccine carrier AOV

Product Code :
LBNY-0001-2100027



Description :

Vaccine storage capacity: 1.5 liters
With lid and shoulder strap with flapper.
Cold life 34 hours at 43°C without openings.
Supplied with 4 x 0.6 liter icepacks as a standard.
Weight fully loaded: 8 kg.

Technical Specification :

Dimensions:
Internal dimensions (LxWxH): 9.40 x 9.40 x 17 cm
External dimensions (LxWxH): 31 x 31 x 30 cm



Laboratory Instrument India