



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Vaccine carrier AOV

**Product Code :**  
LBNY-0001-2100027



**Description :**

Vaccine storage capacity: 1.5 liters  
With lid and shoulder strap with flapper.  
Cold life 34 hours at 43°C without openings.  
Supplied with 4 x 0.6 liter icepacks as a standard.  
Weight fully loaded: 8 kg.

**Technical Specification :**

Dimensions:  
Internal dimensions (LxWxH): 9.40 x 9.40 x 17 cm  
External dimensions (LxWxH): 31 x 31 x 30 cm



**Laboratory Instrument India**