



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Cold box, AOV

**Product Code :**  
LBNY-0001-2100028



**Description :**

Cold box, Large, Long range  
Vaccine storage capacity: 18 litres  
With hinged lid and carrying handle.  
Cold life at 43°C without openings: 126 hours 32 minutes  
Supplied with 50 x 0.3L icepacks as a standard.  
Weight fully loaded: 45.86 kg.

**Technical Specification :**

Dimensions:  
Internal dimensions (LxWxH cm): 53.00 X 37.50 X 30.50  
External dimensions (LxWxH cm): 70.00 X 61.00 X 51.00



**Laboratory Instrument India**