



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cold box, AOV

Product Code :
LBNY-0001-2100028



Description :

Cold box, Large, Long range
Vaccine storage capacity: 18 litres
With hinged lid and carrying handle.
Cold life at 43°C without openings: 126 hours 32 minutes
Supplied with 50 x 0.3L icepacks as a standard.
Weight fully loaded: 45.86 kg.

Technical Specification :

Dimensions:
Internal dimensions (LxWxH cm): 53.00 X 37.50 X 30.50
External dimensions (LxWxH cm): 70.00 X 61.00 X 51.00



Laboratory Instrument India