



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cold box 15 litres

Product Code :
LBNY-0001-2100029



Description :

Vaccine storage capacity: 15 litres
Cold box, Small, Short range
With hinged lid and carrying handles.
Cold life at 43°C without openings: 53 hours 36 minutes
Supplied with 32 x 0.4L icepacks as a standard.
Weight fully loaded: 33.85 Kg

Technical Specification :

Dimensions:
Internal dimensions (LxWxH cm): 48 x 48 x 20
External dimensions (LxWxH cm): 65 x 65 x 37



Laboratory Instrument India