



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Cold box 15 litres

**Product Code :**  
LBNY-0001-2100029



**Description :**

Vaccine storage capacity: 15 litres  
Cold box, Small, Short range  
With hinged lid and carrying handles.  
Cold life at 43°C without openings: 53 hours 36 minutes  
Supplied with 32 x 0.4L icepacks as a standard.  
Weight fully loaded: 33.85 Kg

**Technical Specification :**

Dimensions:  
Internal dimensions (LxWxH cm): 48 x 48 x 20  
External dimensions (LxWxH cm): 65 x 65 x 37



**Laboratory Instrument India**