



Email : sales@laboratoryinstrumentindia.com

Product Name :
Stereo Microscope

Product Code :
LBNY-0003-830005



Description :

Stereo Microscope

Technical Specification :

Stereo Microscope

20X and 40X stereo magnification settings, high quality widefield optical glass lenses, forward widefield binocular head, sturdy pillar stand with metal framework, built-in incident (top) and transmitted (bottom) lights.

Stereo microscope for home or classroom use enables students and hobbyists to inspect biological specimens, rocks, stamps, jewelry, and large specimens that require handling or manipulation.

Binocular viewing head with interchangeable pairs of 10x and 20x widefield eyepieces, adjustable interpupillary distance, and 45-degree inclination to simplify viewing for young users.

Interchangeable 2x and 4x objectives provide low magnification and longer focal length for inspecting large-scale specimens.

Upper and lower illumination has 10W halogen light source to illuminate specimens

Frosted stage plate allows light to pass through from below with transparent specimens, reversible black and white stage plate provides contrast with light- and dark-colored specimens, and stage clips secure the specimen during viewing.

20X-40X widefield stereo magnification settings.

Built-in incident (top) and transmitted (bottom) lights.

Sturdy pillar stand with metal framework.

Backward widefield binocular head.



Laboratory Instrument India