



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cold Box Leff

Product Code :
LBNY-0001-2100030



Description :

Vaccine storage capacity: 15.40 liters
With lid and carrying handle.
Cold life 106 hours 9 minutes at 43°C without openings.
Supplied with 21 x 0.6 liter.
Weight fully loaded: 49.90 kg.

Technical Specification :

Dimensions:
Internal dimensions (LxWxH): 41.5 x 18.5 x 20.0cm
External dimensions (LxWxH): 77 x 54 x 47cm.



Laboratory Instrument India