



Email : sales@laboratoryinstrumentindia.com

Product Name :
VC Nilkamal

Product Code :
LBNY-0001-2100031



Description :

Large, User Independent Freeze Free Vaccine carrier, long range from Nilkamal, for transport of vaccines.
Vaccine storage capacity: 1.51 liters
With lid and shoulder strap.
Cold life at 43°C without openings: 33 hrs 50 min
Supplied with 4 x 0.6 liter.
Weight fully loaded: 7.6 kg.

Technical Specification :

Dimensions:
Internal dimensions (LxWxH): 9.50 x 9.50 x 16.80 cm
External dimensions (LxWxH): 31.8 x 31.8 x 29.5 cm



Laboratory Instrument India