



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
VC Nilkamal

**Product Code :**  
LBNY-0001-2100031



**Description :**

Large, User Independent Freeze Free Vaccine carrier, long range from Nilkamal, for transport of vaccines.  
Vaccine storage capacity: 1.51 liters  
With lid and shoulder strap.  
Cold life at 43°C without openings: 33 hrs 50 min  
Supplied with 4 x 0.6 liter.  
Weight fully loaded: 7.6 kg.

**Technical Specification :**

Dimensions:  
Internal dimensions (LxWxH): 9.50 x 9.50 x 16.80 cm  
External dimensions (LxWxH): 31.8 x 31.8 x 29.5 cm



**Laboratory Instrument India**