



Email : sales@laboratoryinstrumentindia.com

Product Name :
VC BMedical

Product Code :
LBNY-0001-2100032



Description :

Vaccine storage capacity: 1.04 liters
Cold life at 43°C without openings: 32 hrs 07 min
Supplied with 3 x 0.6-liter icepacks (Bmedical WP 0.6L) as a standard.
With lid and shoulder strap
Weight fully loaded: 7.6 kg.

Technical Specification :

Dimensions:
Internal dimensions (LxWxH): 22.3 x 13.9 x 20.8 cm
External dimensions (LxWxH): 34.7 x 28.1 x 43 cm



Laboratory Instrument India