



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Vacc. carr, 1.7 litres

**Product Code :**  
LBNY-0001-2100033



**Description :**

Vaccine storage capacity: 1.7 litres  
Vaccine carrier, Large, Long range  
With lid and shoulder strap.  
Cold life at +43°C without openings: 38 hours 37 minutes  
Supplied with 4 x 0.4L icepacks as a standard.  
Weight fully loaded: 3.34 Kg.

**Technical Specification :**

Dimensions:  
Internal dimensions (LxWxH cm): 17.0 x 17.0 x 16.4  
External dimensions (LxWxH cm): 25.2 x 26.0 x 30.5.



**Laboratory Instrument India**