



Email : sales@laboratoryinstrumentindia.com

Product Name :
Vacc. carr, 1.7 litres

Product Code :
LBNY-0001-2100033



Description :

Vaccine storage capacity: 1.7 litres
Vaccine carrier, Large, Long range
With lid and shoulder strap.
Cold life at +43°C without openings: 38 hours 37 minutes
Supplied with 4 x 0.4L icepacks as a standard.
Weight fully loaded: 3.34 Kg.

Technical Specification :

Dimensions:
Internal dimensions (LxWxH cm): 17.0 x 17.0 x 16.4
External dimensions (LxWxH cm): 25.2 x 26.0 x 30.5.



Laboratory Instrument India