




Email : sales@laboratoryinstrumentindia.com

Product Name :
Vacc.carrier,Gio Style

Product Code :
LBNY-0001-2100034



Description :

Vaccine carrier, Extra-large, Long range
Vaccine storage capacity: 2.6 litres
With lid and carrying handle.
Cold life at +43°C without openings: 32 hours
Supplied with 8 x 0.4L icepacks as a standard.
Weight fully loaded: 6.5 Kg.

Technical Specification :

Dimensions:
Internal dimensions (LxWxH cm): 22 x 18 x 24
External dimensions (LxWxH cm): 29 x 24 x 32



Laboratory Instrument India