




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Cold Box Leff 20 liters

**Product Code :**  
LBNY-0001-2100035



**Description :**

Vaccine storage capacity: 20 liters.  
With lid and two hinged handles.  
Cold life 72 hours 52 minutes at 43°C without openings.  
Supplied with 17 x 0.6 liter icepacks as a standard.  
Weight fully loaded: 49.20 kg.

**Technical Specification :**

Dimensions:  
Internal dimensions (LxWxH): 55.5 x 18.5 x 19.5 cm  
External dimensions (LxWxH): 77.5 x 54.5 x 47.3 cm



**Laboratory Instrument India**