




Email : sales@laboratoryinstrumentindia.com

Product Name :
Cold Box Leff 20 liters

Product Code :
LBNY-0001-2100035



Description :

Vaccine storage capacity: 20 liters.
With lid and two hinged handles.
Cold life 72 hours 52 minutes at 43°C without openings.
Supplied with 17 x 0.6 liter icepacks as a standard.
Weight fully loaded: 49.20 kg.

Technical Specification :

Dimensions:
Internal dimensions (LxWxH): 55.5 x 18.5 x 19.5 cm
External dimensions (LxWxH): 77.5 x 54.5 x 47.3 cm



Laboratory Instrument India