




Email : sales@laboratoryinstrumentindia.com

Product Name :
VC 1.6 litres

Product Code :
LBNY-0001-2100036



Description :

Vaccine storage capacity: 1.6 liters.
With lid and carrying handle.
Cold life 30 hours at 43°C without openings.
Supplied with 4 x 0.4 liter icepacks as a standard.
Weight fully loaded: 6.4 kg.

Technical Specification :

Dimensions:
Internal dimensions (LxWxH): 10 x 10 x 16 cm
External dimensions (LxWxH): 28.5 x 28.5 x 33.5 cm



Laboratory Instrument India