



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Foster WIC Mono Block 10m3

**Product Code :**  
LBNY-0002-210001



**Description :**

Foster Walk-in type cold room, mono block, 10m3 for storage of vaccines at national and sub-national levels.  
Internal dimensions: L=2100 x W=2100 x H=2100 mm  
External dimensions: L=2300 x W=2300 x H=2300 mm

**Technical Specification :**

The electrical control is fitted with phase failure and voltage spike protection.  
Prefabricated Cold room, walk-in type, mono block, supplied with all necessary mechanical and electrical fittings and equipment for quick and easy installation in the field.  
The complete unit is also designed for installation in a housed area such as a warehouse in order to meet temperature design standards.  
The unit is designed to operate from a mains supply of 220/240 V/1/50Hz with a dual refrigeration unit, each refrigeration unit is sized to give 100% stand-by capacity designed to operate at +2°C - +8°C in a maximum ambient of 43°C.  
Units operating from a mains supply of 220V/1 phase/60Hz is also available as an option.



**Laboratory Instrument India**

