



Email : sales@laboratoryinstrumentindia.com

Product Name :
Foster WIC Mono Block 10m3

Product Code :
LBNY-0002-210001



Description :

Foster Walk-in type cold room, mono block, 10m3 for storage of vaccines at national and sub-national levels.
Internal dimensions: L=2100 x W=2100 x H=2100 mm
External dimensions: L=2300 x W=2300 x H=2300 mm

Technical Specification :

The electrical control is fitted with phase failure and voltage spike protection.
Prefabricated Cold room, walk-in type, mono block, supplied with all necessary mechanical and electrical fittings and equipment for quick and easy installation in the field.
The complete unit is also designed for installation in a housed area such as a warehouse in order to meet temperature design standards.
The unit is designed to operate from a mains supply of 220/240 V/1/50Hz with a dual refrigeration unit, each refrigeration unit is sized to give 100% stand-by capacity designed to operate at +2°C - +8°C in a maximum ambient of 43°C.
Units operating from a mains supply of 220V/1 phase/60Hz is also available as an option.



Laboratory Instrument India

