



Email : sales@laboratoryinstrumentindia.com

Product Name :
Foster WIC Mono Block 30m3

Product Code :
LBNY-0002-210002



Description :

Foster Walk-in type cold room, mono block, 30m3 for storage of vaccines at national and sub-national levels.
Internal dimensions: L=3690 x W=3690 x H=2100 mm
External dimensions: L=3890 x W=3890 x H=2300 mm

Technical Specification :

Prefabricated Cold room, walk-in type, mono block, supplied with all necessary mechanical and electrical fittings and equipment for quick and easy installation in the field.
The complete unit is also designed for installation in a housed area such as a warehouse in order to meet temperature design standards.
The unit is designed to operate from a mains supply of 220/240 V/3/50Hz with a dual refrigeration unit, each refrigeration unit is sized to give 100% stand-by capacity designed to operate at +2°C - +8°C in a maximum ambient of 43°C.
Units operating from a mains supply of 220V/3 phase/60Hz is also available as an option.
The electrical control is fitted with phase failure and voltage spike protection.



Laboratory Instrument India

