



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Foster WIC Mono Block 40m3

**Product Code :**  
LBNY-0002-210003



**Description :**

Foster Walk-in type cold room, mono block, 40m3 for storage of vaccines at national and sub-national levels.  
Internal dimensions: L=4260 x W=4260 x H=2100 mm  
External dimensions: L=4460 x W=4460 x H=2300 mm

**Technical Specification :**

The unit is designed to operate from a mains supply of 220/240 V/3/50Hz with a dual refrigeration unit, each refrigeration unit is sized to give 100% stand-by capacity designed to operate at +2°C - +8°C in a maximum ambient of 43°C.

Units operating from a mains supply of 220V/3 phase/60Hz is also available as an option.

Prefabricated Cold room, walk-in type, mono block, supplied with all necessary mechanical and electrical fittings and equipment for quick and easy installation in the field.

The complete unit is also designed for installation in a housed area such as a warehouse in order to meet temperature design standards.

The electrical control is fitted with phase failure and voltage spike protection.



**Laboratory Instrument India**

