



Email : sales@laboratoryinstrumentindia.com

Product Name :
Icepack, 0.3-liter capacity

Product Code :
LBNY-0003-210001



0.3 litres

Description :

Capacity: 0.3 liters
Icepack with screw cap (Type2)
Water content: 0.25 to 0.3 liters

Technical Specification :

External dimensions: 163 x 90 x 34 mm, +/- 2 mm (Type 2)
Removable cap, with internal water seal to prevent leakage.
Filling line indicated on one side.
Walls to be reinforced to prevent swelling.
Supplied empty.



Laboratory Instrument India