



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Icepack, 0.4-liter capacity

**Product Code :**  
LBNY-0003-210002



**0.4 litres**

**Description :**

Icepack with screw cap  
Water content: 0.35 to 0.4 liters  
Capacity: 0.4 liters  
External dimensions: 163 x 94 x 34 mm, +/- 2 mm  
Removable cap, with internal water seal to prevent leakage.  
Filling line indicated on one side.  
Walls to be reinforced to prevent swelling  
Supplied empty.

**Technical Specification :**

Icepack, 0.4-liter capacity



**Laboratory Instrument India**