



Email : sales@laboratoryinstrumentindia.com

Product Name :
Icepack, 0.4-liter capacity

Product Code :
LBNY-0003-210002



0.4 litres

Description :

Icepack with screw cap
Water content: 0.35 to 0.4 liters
Capacity: 0.4 liters
External dimensions: 163 x 94 x 34 mm, +/- 2 mm
Removable cap, with internal water seal to prevent leakage.
Filling line indicated on one side.
Walls to be reinforced to prevent swelling
Supplied empty.

Technical Specification :

Icepack, 0.4-liter capacity



Laboratory Instrument India