



Email : sales@laboratoryinstrumentindia.com

Product Name :
Icepack, 0.6-liter capacity

Product Code :
LBNY-0003-210003



0.6 litres

Description :

Icepack with screw cap
Capacity: 0.6 litres
External dimensions: 190 x 120 x 34mm, +/- 2mm
Water content: 0.55 to 0.6 litres
Removable cap, with internal water seal to prevent leakage.
Filling line indicated on one side.
Walls to be reinforced to prevent swelling.
Supplied empty.

Technical Specification :

Icepack, 0.6-liter capacity



Laboratory Instrument India