



Email : sales@laboratoryinstrumentindia.com

Product Name :

Water Pack B Medical WP 0.3L

Product Code :

LBNY-0003-210004



Description :

Icepack with screw cap (Type1)
Capacity: 0.3 litres
External dimensions: 173 x 120 x 26mm, +/- 2mm (Type1)
Water content: 0.25 to 0.35 litres
Removable cap, with internal water seal to prevent leakage.
Filling line indicated on one side.
Walls to be reinforced in order to prevent swelling.
Supplied empty.

Technical Specification :

Water Pack B Medical WP 0.3L



Laboratory Instrument India