



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Water Pack B Medical WP 0.3L

**Product Code :**

LBNY-0003-210004



**Description :**

Icepack with screw cap (Type1)  
Capacity: 0.3 litres  
External dimensions: 173 x 120 x 26mm, +/- 2mm (Type1)  
Water content: 0.25 to 0.35 litres  
Removable cap, with internal water seal to prevent leakage.  
Filling line indicated on one side.  
Walls to be reinforced in order to prevent swelling.  
Supplied empty.

**Technical Specification :**

Water Pack B Medical WP 0.3L



**Laboratory Instrument India**