



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Water Pack B Medical WP 0.6L

**Product Code :**  
LBNY-0003-210005



**Description :**

Icepack with screw cap  
Capacity: 0.6 litres  
External dimensions: 190 x 120 x 34mm, +/- 2mm  
Water content: 0.55 to 0.6 litres  
Removable cap, with internal water seal to prevent leakage.  
Filling line indicated on one side.  
Walls to be reinforced to prevent swelling.  
Supplied empty.

**Technical Specification :**

Water Pack B Medical WP 0.6L



**Laboratory Instrument India**