



Email : sales@laboratoryinstrumentindia.com

Product Name :
Water Pack B Medical WP 0.6L

Product Code :
LBNY-0003-210005



Description :

Icepack with screw cap
Capacity: 0.6 litres
External dimensions: 190 x 120 x 34mm, +/- 2mm
Water content: 0.55 to 0.6 litres
Removable cap, with internal water seal to prevent leakage.
Filling line indicated on one side.
Walls to be reinforced to prevent swelling.
Supplied empty.

Technical Specification :

Water Pack B Medical WP 0.6L



Laboratory Instrument India