



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Mains Powered vaccine/water pack freezer 300 L

**Product Code :**  
LBNY-0004-210002



**Description :**

Gros volume-Freezer: 300 L

Volume category: >150L

Equipment operating voltage and frequency: 220-240 Volt 50Hz /60Hz, or 100-127 Volt 50Hz /60Hz, 1ph.

Energy consumption, stable running (kWh/24 hours): 3.54

Operating temperature: Max rated ambient temperature: +43°C

Refrigerant: R600a

Prequalified Mains Powered Vaccine/Water-Pack Freezer for storage of vaccines at -15°C to -25°C and freezing of water packs.

**Technical Specification :**

Mains Powered vaccine/water pack freezer 300 L



**Laboratory Instrument India**