



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mains Powered vaccine/water pack freezer 281 L

Product Code :
LBNY-0004-210004



Description :

Gross volume-Freezer: 281 L
Volume category: >150L
Energy consumption, stable running (kWh/24 hours): 4.23
Operating temperature: Max rated ambient temperature: +43°C
Refrigerant: R600a, 75 grams
Equipment operating voltage and frequency: 220-240 Volt 50Hz /60Hz, or 100-127 Volt 50Hz /60Hz, 1ph.
Prequalified Mains Powered Vaccine/Water-Pack Freezer for storage of vaccines at -15°C to -25°C and freezing of water packs.

Technical Specification :

Mains Powered vaccine/water pack freezer 281 L



Laboratory Instrument India