



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Blood Glucose Meter, with Test Strips

**Product Code :**  
LBNY-0003-830007



**Description :**

Blood Glucose Meter, with Test Strips

**Technical Specification :**

Blood glucose meter, with test strips  
Glucose Test Range - 20-600 mg/dL  
Result Calibration – Plasma  
Test Time- 5 seconds  
Sample Size- 0.6 L  
Operating Temperature - 5°C-45°C  
Operating Humidity- 10-90% RH  
Memory capacity- 200 tests  
Battery Type - 3V Li-Battery X 1  
Battery Life –1,000 tests  
Meter Dimensions – 85 X 52 X 15 mm  
LCD Size - 42 X 39 mm  
Meter Weight - 50g.  
Test strips , pack of 50

---

# Laboratory Instrument India

