



Email : sales@laboratoryinstrumentindia.com

Product Name :
Blood Glucose Meter, with Test Strips

Product Code :
LBNY-0003-830007



Description :

Blood Glucose Meter, with Test Strips

Technical Specification :

Blood glucose meter, with test strips
Glucose Test Range - 20-600 mg/dL
Result Calibration – Plasma
Test Time- 5 seconds
Sample Size- 0.6 L
Operating Temperature - 5°C-45°C
Operating Humidity- 10-90% RH
Memory capacity- 200 tests
Battery Type - 3V Li-Battery X 1
Battery Life –1,000 tests
Meter Dimensions – 85 X 52 X 15 mm
LCD Size - 42 X 39 mm
Meter Weight - 50g.
Test strips , pack of 50

Laboratory Instrument India

