



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mains Powered Vaccine Refrigerator 225 L

Product Code :
LBNY-0004-210006



Description :

Vaccine storage capacity: 225 L
Volume category: >120L
Refrigerant: R600a, 68 grams
Energy consumption, stable running (kWh/24 hours): 2.04
Prequalified Intermittent Mains Powered Refrigerator for storage of vaccines at +2 °C to +8°C.
Operating temperature: Min and Max rated ambient temperature: +10 C and +43°C
UIFPF- Grade A
Equipment operating voltage and frequency: 220-240 Volt, 50Hz/60Hz or 100-127 Volt 50Hz/60Hz, 1ph.

Technical Specification :

Mains Powered Vaccine Refrigerator 225 L



Laboratory Instrument India