



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Mains Powered Combined Vaccine Refrigerator and Water pack Freezer

**Product Code :**  
LBNY-0004-210009



**Description :**

Volume category: Vaccine storage capacity: 52.5 L

Refrigerant: R600a, 60 grams

Energy consumption, stable running (kWh/24 hours): 0.63

Operating temperature: Min and Max rated ambient temperature: +10 °C and +43°C

UIFPF- Grade A

Prequalified Intermittent Mains Powered Combined Refrigerator and Water-Pack Freezer for storage of vaccines at +2 °C to +8°C and freezing of water packs.

220-240 Volt 50Hz /60Hz, or 100-127 Volt 50Hz /60Hz, 1ph.

**Technical Specification :**

Mains Powered Combined Vaccine Refrigerator and Water pack Freezer



**Laboratory Instrument India**