



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mains Powered Vaccine Refrigerator 240 L

Product Code :
LBNY-0004-2100011



Description :

Vaccine storage capacity: 240 L
Volume category: ≥ 120 L
Refrigerant: R600a, 60 grams
Energy consumption, stable running (kWh/24 hours): 0.85.
Equipment operating voltage and frequency: 220-240 Volt 50Hz /60Hz or 115 Volt /60Hz, 1ph.
Operating temperature: Min and Max rated ambient temperature: +5 °C and +43°C
UIFPF- Grade A
PQS Prequalified Intermittent Mains Powered Refrigerator for storage of vaccines at +2 °C to +8°C.

Technical Specification :

Mains Powered Vaccine Refrigerator 240 L



Laboratory Instrument India