



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Mains Powered Vaccine Refrigerator 240 L

**Product Code :**  
LBNY-0004-2100011



**Description :**

Vaccine storage capacity: 240 L  
Volume category:  $\geq 120$  L  
Refrigerant: R600a, 60 grams  
Energy consumption, stable running (kWh/24 hours): 0.85.  
Equipment operating voltage and frequency: 220-240 Volt 50Hz /60Hz or 115 Volt /60Hz, 1ph.  
Operating temperature: Min and Max rated ambient temperature: +5 °C and +43°C  
UIFPF- Grade A  
PQS Prequalified Intermittent Mains Powered Refrigerator for storage of vaccines at +2 °C to +8°C.

**Technical Specification :**

Mains Powered Vaccine Refrigerator 240 L



**Laboratory Instrument India**