



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mains Powered Vaccine Refrigerator 184 L

Product Code :
LBNY-0004-2100012



Description :

Vaccine storage capacity: 184 L
Volume category: $\geq 120L$
Refrigerant: R600a, 70 grams
Energy consumption, stable running (kWh/24 hours): 0.803
Operating temperature: Min and Max rated ambient temperature: +5 °C and +43°C
UIFPF- Grade A
PQS Prequalified Intermittent Mains Powered Refrigerator for storage of vaccines at +2 °C to +8°C.
Equipment operating voltage and frequency.
220-240 Volt, 50Hz/60Hz or 100-127 Volt 50Hz/60Hz, 1ph.

Technical Specification :

Mains Powered Vaccine Refrigerator 184 L



Laboratory Instrument India