



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Mains Powered Vaccine Refrigerator 61 L

**Product Code :**  
LBNY-0004-2100013



**Description :**

Vaccine storage capacity: 61L  
Volume category:  $\geq 60$ , Refrigerant: R600a, 60 grams  
Equipment operating voltage and frequency.  
220-240 Volt and frequency 50Hz/60Hz, 1ph.  
Operating temperature: Min and Max rated ambient temperature: +5 °C and +43°C  
UIFPF- Grade A  
Prequalified Intermittent Mains Powered Refrigerator for storage of vaccines at +2°C to +8°C.

**Technical Specification :**

Mains Powered Vaccine Refrigerator 61 L



**Laboratory Instrument India**