



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mains Powered Vaccine Refrigerator 61 L

Product Code :
LBNY-0004-2100013



Description :

Vaccine storage capacity: 61L
Volume category: ≥ 60 , Refrigerant: R600a, 60 grams
Equipment operating voltage and frequency.
220-240 Volt and frequency 50Hz/60Hz, 1ph.
Operating temperature: Min and Max rated ambient temperature: +5 °C and +43°C
UIFPF- Grade A
Prequalified Intermittent Mains Powered Refrigerator for storage of vaccines at +2°C to +8°C.

Technical Specification :

Mains Powered Vaccine Refrigerator 61 L



Laboratory Instrument India