



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mains Powered Vaccine Refrigerator 211 L

Product Code :
LBNY-0004-2100014



Description :

Vaccine storage capacity: 211L
Volume category: $\geq 120L$
Refrigerant: R600a, 80 grams
Operating temperature: Min and Max rated ambient temperature: $+5^{\circ}C$ and $+43^{\circ}C$
UIFPF- Grade A
Prequalified Intermittent Mains Powered Refrigerator for storage of vaccines at $+2^{\circ}C$ to $+8^{\circ}C$.
Equipment operating voltage and frequency: 220-240 Volt and frequency 50Hz / 60Hz, 1ph.

Technical Specification :

Mains Powered Vaccine Refrigerator 211 L



Laboratory Instrument India